

**Homogenisers, Safe-Seal, borosilicate glass 3.3/PTFE****GPE****Homogenisers, Safe-Seal, borosilicate glass 3.3/PTFE**

Supplied complete with a PTFE pestle, stainless steel handle and borosilicate glass mortar which is stoppered in PTFE to protect against splash back during the grinding process. Autoclavable.



Code	Description	Packaging
LLG09164504	Homogenisers, Safe-Seal, borosilicate glass 3.3/PTFE, 30 ml, Pestle L x o.d. 250 x 19 mm, Mortar L x o.d. 180 x 23 mm	1 pz
LLG09164502	Homogenisers, Safe-Seal, borosilicate glass 3.3/PTFE, 10 ml, Pestle L x o.d. 250 x 13 mm, Mortar L x o.d. 150 x 17 mm	1 pz
LLG09164503	Homogenisers, Safe-Seal, borosilicate glass 3.3/PTFE, 15 ml, Pestle L x o.d. 250 x 15 mm, Mortar L x o.d. 155 x 19 mm	1 pz
LLG09164501	Homogenisers, Safe-Seal, borosilicate glass 3.3/PTFE, 5 ml, Pestle L x o.d. 250 x 10 mm, Mortar L x o.d. 140 x 15 mm	1 pz



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

