

**Homogenisers, Griffiths tube, Borosilicate glass 3.3****GPE****Homogenisers, Griffiths tube, Borosilicate glass 3.3**

For grinding bacterial emulsions and tissues. Rounded pestle and mortar ends, ensure perfect grinding. Autoclavable. Large Griffiths tube is ideal for coarse material.



Code	Description	Packaging
LLG09164550	Homogenisers, Griffiths tube, Borosilicate glass 3.3, 5 ml, Pestle L x o.d. 150 x 12 mm, Mortar L x o.d. 95 x 18 mm, small	1 pz
LLG09164551	Homogenisers, Griffiths tube, Borosilicate glass 3.3, 15 ml, Pestle L x o.d. 305 x 14 mm, Mortar L x o.d. 220 x 18 mm, large	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

