


**Homogenisers, Round Body, glass****GPE****Homogenisers, Round Body, glass**

The all-glass Uni-Form range of round body tissue grinders is made from sterilisable heat-resisting glass in capacities 5ml to 55ml. Models are available with rod ends for motor drive, or hollow pestles designed for use with coolants/dry ice. Bodies are available in clear precision glass or ground glass for more effective homogenisation of tougher tissues.

Grinding clearance pestle and mortar 0.152mm to 0.254mm.



Code	Description	Packaging
LLG09164661	Homogenisers, Round Body, glass, 10 ml, Pestle L x o.d. 225 x 13 mm, Mortar L x o.d. 130 x 17 mm, Mortar clear glass, Pestle hollow	1 pz
LLG09164666	Homogenisers, Round Body, glass, 10 ml, Pestle L x o.d. 225 x 13 mm, Mortar L x o.d. 130 x 17 mm, Mortar clear glass, Pestle with rod end	1 pz
LLG09164671	Homogenisers, Round Body, glass, 10 ml, Pestle L x o.d. 225 x 13 mm, Mortar L x o.d. 130 x 17 mm, Mortar ground glass, Pestle hollow	1 pz
LLG09164676	Homogenisers, Round Body, glass, 10 ml, Pestle L x o.d. 225 x 13 mm, Mortar L x o.d. 130 x 17 mm, Mortar ground glass, Pestle with rod end	1 pz
LLG09164662	Homogenisers, Round Body, glass, 15 ml, Pestle L x o.d. 240 x 15 mm, Mortar L x o.d. 140 x 19 mm, Mortar clear glass, Pestle hollow	1 pz
LLG09164667	Homogenisers, Round Body, glass, 15 ml, Pestle L x o.d. 240 x 15 mm, Mortar L x o.d. 140 x 19 mm, Mortar clear glass, Pestle with rod end	1 pz
LLG09164672	Homogenisers, Round Body, glass, 15 ml, Pestle L x o.d. 240 x 15 mm, Mortar L x o.d. 140 x 19 mm, Mortar ground glass, Pestle hollow	1 pz
LLG09164677	Homogenisers, Round Body, glass, 15 ml, Pestle L x o.d. 240 x 15 mm, Mortar L x o.d. 140 x 19 mm, Mortar ground glass, Pestle with rod end	1 pz
LLG09164663	Homogenisers, Round Body, glass, 30 ml, Pestle L x o.d. 295 x 19 mm, Mortar L x o.d. 170 x 23 mm, Mortar clear glass, Pestle hollow	1 pz
LLG09164668	Homogenisers, Round Body, glass, 30 ml, Pestle L x o.d. 295 x 19 mm, Mortar L x o.d. 170 x 23 mm, Mortar clear glass, Pestle with rod end	1 pz
LLG09164673	Homogenisers, Round Body, glass, 30 ml, Pestle L x o.d. 295 x 19 mm, Mortar L x o.d. 170 x 23 mm, Mortar ground glass, Pestle hollow	1 pz
LLG09164678	Homogenisers, Round Body, glass, 30 ml, Pestle L x o.d. 295 x 19 mm, Mortar L x o.d. 170 x 23 mm, Mortar ground glass, Pestle with rod end	1 pz
LLG09164660	Homogenisers, Round Body, glass, 5 ml, Pestle L x o.d. 205 x 10 mm, Mortar L x o.d. 105 x 15 mm, Mortar clear glass, Pestle hollow	1 pz
LLG09164665	Homogenisers, Round Body, glass, 5 ml, Pestle L x o.d. 205 x 10 mm, Mortar L x o.d. 105 x 15 mm, Mortar clear glass, Pestle with rod end	1 pz



Code	Description	Packaging
LLG09164670	Homogenisers, Round Body, glass, 5 ml, Pestle L x o.d. 205 x 10 mm, Mortar L x o.d. 105 x 15 mm, Mortar ground glass, Pestle hollow	1 pz
LLG09164675	Homogenisers, Round Body, glass, 5 ml, Pestle L x o.d. 205 x 10 mm, Mortar L x o.d. 105 x 15 mm, Mortar ground glass, Pestle with rod end	1 pz
LLG09164664	Homogenisers, Round Body, glass, 55 ml, Pestle L x o.d. 330 x 25 mm, Mortar L x o.d. 205 x 31 mm, Mortar clear glass, Pestle hollow	1 pz
LLG09164669	Homogenisers, Round Body, glass, 55 ml, Pestle L x o.d. 330 x 25 mm, Mortar L x o.d. 205 x 31 mm, Mortar clear glass, Pestle with rod end	1 pz
LLG09164674	Homogenisers, Round Body, glass, 55 ml, Pestle L x o.d. 330 x 25 mm, Mortar L x o.d. 205 x 31 mm, Mortar ground glass, Pestle hollow	1 pz
LLG09164679	Homogenisers, Round Body, glass, 55 ml, Pestle L x o.d. 330 x 25 mm, Mortar L x o.d. 205 x 31 mm, Mortar ground glass, Pestle with rod end	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

