

**Filter Units Nalgene, Low Profile, PES Membrane, sterile****Thermo Scientific****Filter Units Nalgene, Low Profile, PES Membrane, sterile**

Compact Nalgene filter units with PES membrane are ideally suited for quick filtration of samples between 50 to 100 ml.

- Unitary construction with non-separable upper and receiver components
- Padless membrane support minimizes foaming of proteinaceous samples
- Pyrogen-free
- Sterile



Code	Description	Packaging
LLG06234424	Filter Units, Low Profile, PES Membrane, sterile, Type 524, Volume upper/lower part 115/150 ml, Membrane diam. 50 mm, Pore size 0.2 µm	1 pz.
LLG09046020	Filter Units, Low Profile, PES Membrane, sterile, Type 124, Volume upper/lower part 115/150 ml, Membrane diam. 50 mm, Pore size 0.45 µm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

