

**Photometers AL100****Lovibond****Photometers AL100**

The AQUALYTIC® AL100 Photometers are designed for analysis of different parameters in drinking water, industrial process water and waste water. With a total of more than 30 different types, instruments are available for routine testing as well as for special applications. The AL100 uses high quality interference filters with long-life LEDs as a light source without any moving parts in a transparent sample chamber. The units supply accurate, reproducible results very quickly. Other major advantages include ease of operation, ergonomic design, compact dimensions and safe handling. Provided with calibration and software-based adjustment options the AL100 is suitable for use as a testing instrument. The tests are conducted using either AQUALYTIC® tablet reagents with long-term stability and a shelf life of 5 or 10 years, VARIO powder reagents or using liquid reagents.

- Optics: High-quality Interference Filters
- Wavelength accuracy: ± 1 nm.
- Power supply: 4 micro batteries (AAA), capacity approx. 17 hours or 5000 tests.
- Automatic switch-off.
- LCD Display (backlit on keypress).
- Internal ring memory for 16 data sets.
- Dimensions (L x W x H): 155 x 75 x 35 mm.



Scope of supply: Photometer in a plastic carrying case, with reagents for 100 tests, cells, batteries and standard accessories.

Optional: IRiM - infrared interface module for data transfer to computer or printer.



Code	Description	Packaging
LLG07910364	Photometers AL100, 2 in 1 Chloride/pH with liquid reagents	1 pz.
LLG09699288	Photometers AL100, 3 in 1 Chlorine LR / Chlorine HR / Chlorine dioxide tablet reagents	1 pz.
LLG07627352	Reference standard kit Chloride 0,2-1,0 mg	1 pz.
LLG09699280	Photometers AL100, Ammonium tablet reagents	1 pz.
LLG09699291	Photometers AL100, Reference standard kit pH	1 pz.
LLG09699281	Photometers AL100, Chlorine dioxide tablet reagents	1 pz.
LLG06252939	Photometers AL100, Chlorine liquid reagents	1 pz.
LLG06251486	Photometers AL100, Chlorine powder reagents	1 pz.
LLG09699283	Photometers AL100, Copper tablet reagents	1 pz.
LLG09699282	Photometers AL100, Iron tablet reagents	1 pz.
LLG09699284	Photometers AL100, Manganese LR tablet reagents	1 pz.
LLG09699285	Photometers AL100, Phosphate tablet reagents	1 pz.
LLG06252821	Photometers AL100, Reference standard kit Chloride	1 pz.
LLG09699289	Photometers AL100, 5 in 1 Chlorine / pH / Cyanuric acid / Alkalinity-M / Calcium hardness tablet reagents	1 pz.





CARLO ERBA Reagents S.A.S.
Chaussée du Vexin – 27106 Val
de Reuil cedex
N° TVA: FR 63391048824
Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



*CARLO ERBA
Reagents operates with
a Certified Quality
Management System*

