


**Tumbling shaker GyroTwister / GyroTwister Adjustable****Corning BV****Tumbling shaker GyroTwisterä / GyroTwisterä Adjustable**

- Gyratory action provides gentle, thorough mixing
- Ideal for staining gels, washing membranes and mixing samples in tubes
- Fixed angle and adjustable tilt models
- Safe for cold room or incubator use
- Platform: 30cm x 30cm
- Speed range: 3-60rpm
- Timer: 20min to 120min or continuous
- Motion/Pitch: tumbling/adjustable +/- 10°
- Motion/Pitch: tumbling/Fixed 7° pitch
- Max. load: 5kg
- Operating range: +4 to 70 °C
- Dimensions: 33cm x 35.2cm x 21cm
- Weight: 6.6kg
- Electrical: 230V, 50/60Hz

Shaker comes with 30cm x 30cm non-slip platform and dimpled tube mat.



Code	Description	Packaging
LLG06227093	GyroTwister adjustable shaker with 30x30 cm non-slip platform 230 V, incl. EU & UK power cords	1 pz.
LLG09595335	Accessories for Tumbling shaker GyroTwister / GyroTwister Adjustable, Staking platform 30cm x 30cm, with hardware	1 pz.
LLG09595338	Accessories for Tumbling shaker GyroTwister / GyroTwister Adjustable, Elasticated tie-down	4 pz.
LLG09595336	Accessories for Tumbling shaker GyroTwister / GyroTwister Adjustable, Large platform, 40 x 40cm, with non-slip rubber mat	1 pz.
LLG09595337	Accessories for Tumbling shaker GyroTwister / GyroTwister Adjustable, Extra large platform, 50cm x 50cm, with non-slip rubber mat	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

