

**Radio controlled alarm clock****TFA Dostmann****Radio controlled alarm clock**

Analogue display with pointer, digital display of seconds, alarm with snooze function, backlight, battery 1.5 V AA included.

Dimensions (W x H x D): 96 x 116 x 55 mm

Weight: 140g



Code	Description	Packaging
LLG09262540	Radio controlled alarm clock, Type Radio controlled alarm clock, Measuring capacity 23:59:59 hr: min: sec	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

