

**LLG Laboratory Chairs, PU Foam****LLG Labware****LLG Laboratory Chairs, PU Foam**

According to EN 1728, EN 1021-1 and 1021-2, EN 61340-5-1. Surface enables easy cleaning and disinfecting. Revolving, height adjustable seat. Seat depth and angle and height of the backrest can be adjusted which always guarantees an ergonomic seating position. Metal 5-star bar with 5 castors for soft and 5 castors for hard and soft floors or gliders. Armrests and chromed footring as outlined below.



Code	Description	Packaging
LLG09732201	LLG-Laboratory chair PUR Standard IIPU	1 pz
LLG09732200	LLG Laboratory Chairs, PU Foam, Type Stuhl PUR Standard II plus, 5 gliders + footring + armrests	1 pz
LLG09732213	LLG Laboratory Chairs, PU Foam, Type Stop-and-go Rollen, for hard floors	5 pz.
LLG09732203	LLG-Laboratory chair PUR Standard IPU f	1 pz
LLG09732202	LLG-Laboratory chair PUR Standard I plu	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

