


**Beakers, glass, jacketed with PTFE-Olive, type T-GL****KGW****Beakers, glass, jacketed with PTFE-Olive, type T-GL**

DURAN®. For heating and cooling liquids. Thermostating liquid is circulating between the double walls of the beaker jacket. With glass-screw connections GL14 with PTFE olives (O.D. 9mm, T-GL2000 with GL18 and O.D. 11mm) for thermostat hoses. Coolant pressure: 0.5bar at maximum, operating temperature: +120°C at maximum. Special sizes available on request.



Code	Description	Packaging
LLG09032377	Beakers, glass, jacketed wiht PTFE-Olive, type T-GL, 2000 ml, Neck thread 18 GL, Int. diam. 112 mm, Ext. diam. 135 mm, Int. height 200 mm	1 pz.
LLG09032372	Beakers, glass, jacketed wiht PTFE-Olive, type T-GL, 250 ml, Neck thread 14 GL, Int. diam. 55 mm, Ext. diam. 70 mm, Int. height 115 mm	1 pz.
LLG09032374	Beakers, glass, jacketed wiht PTFE-Olive, type T-GL, 600 ml, Neck thread 14 GL, Int. diam. 77 mm, Ext. diam. 95 mm, Int. height 135 mm	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

