

**Dilution microliter pipette Acura manual 810, fix****Socorex****Dilution microliter pipette Acura® manual 810, fix**

This air displacement micropipette allows subsequent pipetting of 1 and 0.1 mL of same liquid.

The two pre-calibrated fixed volumes eliminate any setting. A universal metal nozzle accommodates straw tips (length 190 mm, diameter 4 mm), most adequate to aspirate from deep, narrow vials or Stomacher bags.

Ideal alternative to conventional graduated glass pipettes when performing serial 1:10 dilutions in bacteriology.



Code	Description	Packaging
LLG06203235	Dilution pipette Acura 810 fix, 1/0.1 ml	1 pz



CARLO ERBA Reagents S.A.S.
Chaussée du Vexin – 27106 Val
de Reuil cedex
N° TVA: FR 63391048824
Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
Reagents operates with
a Certified Quality
Management System

