

**LED-Flashlight, Lumatic Spot****TFA Dostmann****LED-Flashlight, Lumatic Spot**

3 Watt CREE LED, range up to 120m, zoom or focus for diffused light. Durable aluminum housing, splash-proof, with carrying strap and bag, batteries 3 x 1.5V AAA.

Dimensions (Dia. x L):

Weight: 34 x 115mm

170g



Code	Description	Packaging
LLG09156233	LED-Flashlight, Lumatic Spot, Type Lumatic Spot	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

