




Upright Freezer Racks

Upright Freezer Racks

Stainless steel racks, suitable for all brands of upright freezers. For boxes up to 136 mm x 136 mm.



Code	Description	Packaging
LLG09698971	Upright Freezer Racks, Type 5 x 4, Width 558 mm, Depth 139 mm, Height 276 mm, Height Box approx. 53 mm	1 pz.
LLG09698975	Upright Freezer Racks, Type 3 x 4, Width 556 mm, Depth 139 mm, Height 240 mm, Height Box approx. 78 mm	1 pz.
LLG09698976	Upright Freezer Racks, Type 3 x 3, Width 419 mm, Depth 139 mm, Height 240 mm, Height Box approx. 78 mm	1 pz.
LLG09698960	Upright Freezer Racks, Type 4 x 4, Width 556 mm, Depth 140 mm, Height 326 mm, Height Box approx. 78 mm	1 pz.
LLG09698961	Upright Freezer Racks, Type 5 x 4, Width 556 mm, Depth 140 mm, Height 407 mm, Height Box approx. 78 mm	1 pz.
LLG09698950	Upright Freezer Racks, Type 4 x 3, Width 140 mm, Depth 419 mm, Height 221 mm, Height Box approx. 50 mm	1 pz
LLG09698954	Upright Freezer Racks, Type 2 x 3, Width 140 mm, Depth 423 mm, Height 164 mm, Height Box approx. 75 mm	1 pz.
LLG09698959	Upright Freezer Racks, Type 6 x 4, Width 140 mm, Depth 562 mm, Height 331 mm, Height Box approx. 50 mm	1 pz.
LLG09698974	Upright Freezer Racks, Type 5 x 3, Width 419 mm, Depth 139 mm, Height 276 mm, Height Box approx. 53 mm	1 pz.
LLG09698970	Upright Freezer Racks, Type 3 x 4, Width 558 mm, Depth 139 mm, Height 168 mm, Height Box approx. 53 mm	1 pz.
LLG09698973	Upright Freezer Racks, Type 4 x 5, Width 688 mm, Depth 139 mm, Height 221 mm, Height Box approx. 53 mm	1 pz.
LLG09698969	Upright Freezer Racks, Type 4 x 4, Width 558 mm, Depth 139 mm, Height 221 mm, Height Box approx. 53 mm	1 pz.
LLG09698972	Upright Freezer Racks, Type 7 x 4, Width 556 mm, Depth 140 mm, Height 392 mm, Height Box approx. 53 mm	1 pz.





CARLO ERBA Reagents S.A.S.
Chaussée du Vexin – 27106 Val
de Reuil cedex
N° TVA: FR 63391048824
Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



*CARLO ERBA
Reagents operates with
a Certified Quality
Management System*

