

**Condensers, ground glass joint, Allihn, borosilicate glass 3.3****Condensers, ground glass joint, Allihn, borosilicate glass 3.3**

Borosilicate glass 3.3 which is resistant to heat and almost all chemicals. They meet the highest quality standards according to DIN 12581. Polypropylene side arm has a screw thread connector for easy and safe fitting of the tubing. The side arms accept tubings with 8/9mm bore size.



Code	Description	Packaging
LLG04008406	Condensers, ground glass joint, Allihn, borosilicate glass 3.3, with PP side arm, Effective length 160 mm, NS 14/23	1 pz.
LLG04008408	Condensers, ground glass joint, Allihn, borosilicate glass 3.3, with PP side arm, Effective length 400 mm, NS 29/32	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

