



Condensers, ground glass joint, Allihn, borosilicate glass 3.3

Condensers, ground glass joint, Allihn, borosilicate glass 3.3

Borosilicate glass 3.3 which is resistant to heat and almost all chemicals. They meet the highest quality standards according to DIN 12581. Polypropylene side arm has a screw thread connector for easy and safe fitting of the tubing. The side arms accept tubings with 8/9mm bore size.



| Code | Description | Packaging |
|-------------|---|-----------|
| LLG04008406 | Condensers, ground glass joint, Allihn, borosilicate glass 3.3, with PP side arm, Effective length 160 mm, NS 14/23 | 1 pz. |
| LLG04008408 | Condensers, ground glass joint, Allihn, borosilicate glass 3.3, with PP side arm, Effective length 400 mm, NS 29/32 | 1 pz. |



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

