

**Retort stand base, tripod pattern****Bochem****Retort stand base, tripod pattern**

Cast iron varnished. Thread M10. Laboratory stands are supplied without support rod. Please select the suitable rod from "support rods".



Code	Description	Packaging
LLG06401256	Retort stand base, tripod pattern, Foot length 105 mm, Weight 1200 g	1 pz.
LLG06084197	Retort stand base, tripod pattern, Foot length 145 mm, Weight 2000 g	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

