

**Staining jar, PMP (TPX), Schiefferdecker pattern****Kartell****Staining jar, PMP (TPX®), Schiefferdecker pattern**

Accepts 20 slides 76mm x 26mm (10 pairs back to back). Transparent. Complete with lid.

Note: Xylene can cause etching after 2 days at room temperature.



Code	Description	Packaging
LLG09161352	Staining jar, PMP (TPX®), Schiefferdecker pattern, Width 65 mm, Length 76 mm, Height 45 mm	1 pz



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

