


**Beaker tongs, chromated steel****ISOLAB****Beaker tongs, chromated steel**

PVC coated Insulated jaws are used to hold all types of beakers.



Code	Description	Packaging
LLG06237287	Beaker tongs, chromated steel, Length 300 mm, For beakers volume up to 1000ml, Temp. range 120 max. °C	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

