


**Pycnometers, Borosilicate glass 3.3., calibrated****ISOLAB****Pycnometers, Borosilicate glass 3.3., calibrated**

ISO 3507. Individually calibrated to contain (In) to 0.001ml. Net volume capacity is engraved on the body. Each bottle is calibrated together with its stopper. In order to differentiate easily, both the body and the stopper have the same code numbers printed on themselves.



Code	Description	Packaging
LLG06241602	Pycnometers, Borosilicate glass 3.3., calibrated, 10 ml	1 pz.
LLG04008329	Pycnometers, Borosilicate glass 3.3., calibrated, 100 ml	1 pz.
LLG04008328	Pycnometers, Borosilicate glass 3.3., calibrated, 25 ml	1 pz.
LLG04008327	Pycnometers, Borosilicate glass 3.3., calibrated, 5 ml	1 pz.
LLG06241603	Pycnometers, Borosilicate glass 3.3., calibrated, 50 ml	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

