

**Tripod stand****ISOLAB****Tripod stand**

Chrome plated steel. Tripods are meant to be used in conjunction with wire gauze, triangles or a sand bath. Holding ring with 150 mm diameter can be reduced to 125 mm with the supplied reduction attachment.



Code	Description	Packaging
LLG04008512	Tripod stand, Diam. 150 mm, Adjustable height 200 ... 300 mm	1 pz.
LLG06242704	Tripod stand, Diam. 125 mm, Total height 200 mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

