

**LLG-Measuring jugs with handle, PP****LLG Labware****LLG-Measuring jugs with handle, PP**

Made of highly transparent Polypropylen. Functional spout to prevent dripping.
With easy-to-read blue graduation for precise filling with liquids.

- Autoclavable at 121 °C
- Open handle
- Stackable
- According to ISO 7056



| Code | Description | Packaging |
|-------------|---|-----------|
| LLG09013680 | Measuring jugs with handle, PP, Nominal capacity 250 ml, Grad. 5 ml | 2 pz |
| LLG09013681 | Measuring jugs with handle, PP, Nominal capacity 500 ml, Grad. 10 ml | 2 pz |
| LLG09013682 | Measuring jugs with handle, PP, Nominal capacity 1000 ml, Grad. 10 ml | 2 pz |
| LLG09013684 | Measuring jugs with handle, PP, Nominal capacity 3000 ml, Grad. 25 ml | 1 pz |
| LLG09013683 | Measuring jugs with handle, PP, Nominal capacity 2000 ml, Grad. 20 ml | 1 pz |



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

