



Stirring heating mantles LabHEAT KM-MER series

SAF Wärmetechnik

Stirring heating mantles LabHEAT® KM-MER series

LabHEAT® - Heating mantles for flasks, round bottom. Designed and technical identical to series KM-ME, but **in addition equipped with a magnetic stirrer**; rotational frequency up to 1600rpm max. heating element temperature 450°C, 1.5m power supply (earthed) cable and RCD (residual current detection); nominal voltage 230V AC.



Code	Description	Packaging
LLG09642545	Stirring heating mantles LabHEAT® KM-MER series, 100 ml, Flask diam. 64 mm, Rating 100 W, Heating zones 1	1 pz.
LLG09642548	Stirring heating mantles LabHEAT® KM-MER series, 1000 ml, Flask diam. 131 mm, Rating 300 W, Heating zones 2	1 pz.
LLG09642546	Stirring heating mantles LabHEAT® KM-MER series, 250 ml, Flask diam. 85 mm, Rating 150 W, Heating zones 2	1 pz.
LLG09642547	Stirring heating mantles LabHEAT® KM-MER series, 500 ml, Flask diam. 105 mm, Rating 200 W, Heating zones 2	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

