

**LLG-Wide-mouth bottles, with screw cap, LDPE****LLG Labware****LLG-Wide-mouth bottles, with screw cap, LDPE**

Natural, translucent. Thick-walled. Ideal for taking, storing and transporting samples of liquids, pastes and granulated products.
Ideal for liquids, pastes and granulates.



| Code | Description | Packaging |
|-------------|---|-----------|
| LLG09073510 | Wide-mouth bottles, with screw cap, LDPE, 50 ml, Diam. 39 mm, Thread external diam. 32 mm, Height 76 mm | 50 pz |
| LLG09073511 | Wide-mouth bottles, with screw cap, LDPE, 100 ml, Diam. 48 mm, Thread external diam. 32 mm, Height 93.5 mm | 20 pz |
| LLG09073512 | Wide-mouth bottles, with screw cap, LDPE, 250 ml, Diam. 64 mm, Thread external diam. 40 mm, Height 120 mm | 10 pz |
| LLG09073513 | Wide-mouth bottles, with screw cap, LDPE, 500 ml, Diam. 77 mm, Thread external diam. 50 mm, Height 154 mm | 10 pz |
| LLG09073514 | Wide-mouth bottles, with screw cap, LDPE, 1000 ml, Diam. 95 mm, Thread external diam. 65 mm, Height 205 mm | 4 pz |
| LLG09073516 | Wide-mouth bottles, with screw cap, LDPE, 2000 ml, Diam. 120 mm, Thread external diam. 65 mm, Height 247 mm | 1 pz |



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

