


**Petri Dishes, PS****Thermo Scientific****Petri Dishes, PS**

Suitable for culturing of fungi, bacteria and other microorganisms. Contact dish (with grid) useful for sampling in hospital environments (for research use only), as well as in the food science and the pharmaceutical industry. All dishes perform well in automatic dispensers due to complete flatness and uniform height.

Material: Polystyrene



Code	Description	Packaging
LLG04008582	Petri Dishes, PS, 250 ml, Working volume 35 ml, Ext. diam. 140 mm, Height 20 mm, Culture surface 145.0 cm ² , Sterile +	



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

