



NS-Glass stoppers, hollow borosilicate glass 3.3

NS-Glass stoppers, hollow borosilicate glass 3.3

Manufactured from borosilicate glass 3.3 which is resistant to heat and almost all chemicals. They meet the highest quality standards according to DIN 12252.

Lower closed end of stopper features a "drip-tip" to ensure adhered liquid is returned to the flask. Hexagonal flat head prevents the stopper from rolling and allows stopper to be stood on end preventing contamination of the ground surface.



Code	Description	Packaging
LLG06258415	NS-Glass stoppers, hollow borosilicate glass 3.3, amber, NS 24/29	1 pz.
LLG06258416	NS-Glass stoppers, hollow borosilicate glass 3.3, amber, NS 12/21	1 pz.
LLG06258417	NS-Glass stoppers, hollow borosilicate glass 3.3, amber, NS 34/35	1 pz.
LLG07629440	NS-Glass stoppers, hollow borosilicate glass 3.3, amber, NS 29/32	1 pz.
LLG07900484	NS-Glass stoppers, hollow borosilicate glass 3.3, amber, NS 10/19	1 pz.
LLG07970573	NS-Glass stoppers, hollow borosilicate glass 3.3, amber, NS 19/26	1 pz.
LLG06240060	NS-Glass stoppers, hollow borosilicate glass 3.3, clear, NS 34/35	1 pz.
LLG06240333	NS-Glass stoppers, hollow borosilicate glass 3.3, clear, NS 24/29	1 pz.
LLG06242396	NS-Glass stoppers, hollow borosilicate glass 3.3, clear, NS 12/21	1 pz.
LLG07900483	NS-Glass stoppers, hollow borosilicate glass 3.3, clear, NS 10/19	1 pz.
LLG09229990	NS-Glass stoppers, hollow borosilicate glass 3.3, clear, NS 14/23	1 pz.
LLG09229991	NS-Glass stoppers, hollow borosilicate glass 3.3, clear, NS 19/26	1 pz.
LLG09229992	NS-Glass stoppers, hollow borosilicate glass 3.3, clear, NS 29/32	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*



