


**C96 MicroWell Plates, PS****Thermo Scientific****C96 MicroWellä Plates, PS**

- C-shaped bottom for optimal washing and liquid retrieval
- FluoroNuncä for fluorescence based immuno assays with minimal auto-fluorescence
- MaxiSorpä and PolySorpä for quantitative and qualitative solid phase immuno assays
- Working volume range: 50 - 250µl/ well



Code	Description	Packaging
LLG09390553	C96 MicroWell Plates, PS, without lid, Colour Clear, Surface MaxiSorp*, Sterile -	
LLG09390554	C96 MicroWell Plates, PS, without lid, Colour Clear, Surface MaxiSorp, Sterile -	
LLG09390557	C96 MicroWell Plates, PS, without lid, Colour Clear, Surface PolySorp*, Sterile -	
LLG06079678	C96 MicroWell Plates, PS, without lid, Colour White, Surface MaxiSorp*, Sterile -	



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

