



Hook connector, steel/stainless steel

Juchheim Laborgeräte GmbH

Hook connector, steel/stainless steel

Hook bosshead steel, epoxy coated or stainless steel with brass pressure block.



Code	Description	Packaging
LLG09223905	Hook connector, steel/stainless steel, M8 stainless steel thumb bolt, Clamping range 13 mm, Material stainless steel	1 pz.
LLG09223906	Hook connector, steel/stainless steel, M8 brass thumb bolt, Clamping range 13 mm, Material steel	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

