




Multi-purpose bosshead light alloy

Juchheim Laborgeräte GmbH

Multi-purpose bosshead light alloy



Code	Description	Packaging
LLG09223904	Multi-purpose bosshead light alloy, with M8 brass thumb bolt, Clamping range 16 mm, Weight 140 g	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

