



Sealing Tapes for MultiWell Plates

Sealing Tapes for MultiWell Plates


Effectively seals all microplate formats for storage, PCR, microscopy, culture and protection. Pre-cut, plate-sized sheets.

Type 1) General use, PCR, storage, photo-sensitive samples

Type 2) Cell Culture, sterile applications

Type 3) Microscopy, good optics, low auto-fluorescence, DMSO compatible

Type 4) Breathable, cell culture



Code	Description	Packaging
LLG09407482	Sealing Tapes for MultiWell Plates, Type 2), Colour Clear, Glue Acrylic, Sterile +, Material Polyester	
LLG06250417	Sealing Tapes for MultiWell Plates, Type 1), Colour Silver, Glue Silicone, Sterile -, Material Aluminium	
LLG04008574	Sealing Tapes for MultiWell Plates, Type 4), Colour White, Glue Acrylate, Sterile +, Material Rayon	



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

