

**Clip-on magnifiers laboCLIP****Eschenbach****Clip-on magnifiers laboCLIP**

The labo-clip system can be easily attached to the spectacle frame. The labo-clip leaves both hands free for working. Comfortable working distance even with a high level of magnification.

- Used with distance spectacles to achieve optimum image quality, fits most eyeglass frames
- When not needed, the lens can be flipped up and out of the way
- Lenses are interchangeable for different magnification needs
- Lens sizes: monocular dia. 40mm, binocular 74.5 x 22.5mm
- Lightweight
- Lens clamp made from glass fibre reinforced plastic
- Supplied with protective box



Code	Description	Packaging
LLG09151171	Clip-on magnifiers laboCLIP, Type labo-clip, binocular, Working Distance 180 mm, Magnification 2.5x/5.0dpt	1 pz
LLG09151175	Clip-on magnifiers laboCLIP, Type labo-clip, binocular, Working Distance 250 mm, Magnification 2x/2.5dpt	1 pz



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

