

**Abbe refractometer AR4****A.KRUSS Optronic****Abbe refractometer AR4**

Classic refractometers for solids and liquids. With integrated digital thermometer. Optionally they can be combined with a thermostat.


- Scales can be adjusted manually using the drive knob
- Measured value can be read via an eyepiece
- Temperature is shown on the LCD display
- Optionally: In combination with the circulation thermostat PT31 the temperature can be maintained at 20 °C or 25 °C

Specifications

Dimensions (W x D x H): 100 x 190 x 270 mm

Weight: 2.5 kg

Power supply: 7.5 V



Code	Description	Packaging
LLG09801110	Abbe refractometer AR4, Type AR4	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

