



Digital hand-held refractometers DR101-60 / DR201-95 / DR-301-95

A.KRÜSS Optronic

Digital hand-held refractometers DR101-60 / DR201-95 / DR-301-95

Digital hand-held refractometer with two different scales, Brix and RI.
Handling is fast and very easy, distilled water is used for calibration. With automatic temperature compensation.

DR301-95 offers 2 additional user-defined scales and can be connected to a computer. Specialities: optional temperature compensation feature and upper and lower tolerance alarm.

Specifications DR101-60 | DR201-95 | DR301-95**Temperature range:** 10 ... 40 °C**Accuracy:** ± 0.5 °C**Temperature compensation:** 10 ... 40 °C**Dimensions:** 110 x 62 x 32 mm | 130 x 80 x 40 mm | 180 x 100 x 60 mm**Weight:** 160 g | 200 g | 500 g**Power supply:** 1.5 V battery | 1.5 V battery | 9 V battery, incl. AC adapter**IP code:** IP65 | IP64 | IP50

Code	Description	Packaging
LLG06236921	Digital hand-held refractometers DR101-60 / DR201-95 / DR-301-95, Type DR301-95	1 pz.
LLG06230350	Digital hand-held refractometers DR101-60 / DR201-95 / DR-301-95, Type DR201-95	1 pz.
LLG09801117	Digital Handrefractometer DR 101-600-60	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

