

**Illuminated pocket magnifiers, mobilux LED****Eschenbach****Illuminated pocket magnifiers, mobilux® LED**

The tried and tested illuminated hand-held magnifiers now feature a new design and improved functions.

- Easy to use: large, easy to operate light switch
- 3 different filters are supplied
- The batteries are easy to change
- Incl. protective case for lens



Code	Description	Packaging
LLG09151181	Illuminated pocket magnifiers, mobilux® LED, dim. Ø 35 mm, Magnification 7x/28.0 dpt	1 pz.
LLG09151182	Illuminated pocket magnifiers, mobilux® LED, dim. Ø 35 mm, Magnification 10x/38.0 dpt	1 pz.
LLG07970900	Illuminated pocket magnifiers, mobilux® LED, dim. Ø 58 mm, Magnification 6x/24.0 dpt	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

