



## Centrifuge tubes, conical tubes, PP/HDPE

### Centrifuge tubes, conical tubes, PP/HDPE

New and improved conical tubes with the largest writing area on the market for increased traceability with additional space for labeling information. Recyclable plastic rack reduces waste in the lab. New conical tubes have higher RCF ratings, which enable more applications from low speed to super speed centrifugation. They meet USP class VI requirements and are non-cytotoxic, pyrogen-free and RNase/DNase free. Graduated, leak proof and sterile (SAL 10<sup>-6</sup>)

Material tubes: PP

Material caps: HDPE



Code	Description	Packaging
LLG09401430	Centrifuge tubes, conical tubes, PP/HDPE, Bulk, 15 ml, Diam. 17 mm, Height 120 mm, Max. 10500 x g* rcf	
LLG09401431	Centrifuge tubes, conical tubes, PP/HDPE, Racked, 15 ml, Diam. 17 mm, Height 120 mm, Max. 10500 x g* rcf	
LLG09401432	Centrifuge tubes, conical tubes, PP/HDPE, Bulk, 50 ml, Diam. 30 mm, Height 115 mm, Max. 17000 x g* rcf	
LLG09401433	Centrifuge tubes, conical tubes, PP/HDPE, Racked, 50 ml, Diam. 30 mm, Height 115 mm, Max. 17000 x g* rcf	



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

