



Heating bath HB 1500 / HB 1500-S

Heju

Heating bath HB 1500 / HB 1500-S

The heating bath HB 1500 is developed for varied tasks. As bath filling water, oils and sand is used. Different sizes are deliverable! The device has an integrated temperature controller with microprocessor. It is freely programmable and has also self optimization function and rampfunction!

- Microprocessor controller with 4-digits 7-segment display, 10 mm high, green,
- With selfoptimization and rampfunction
- Measurement and control range +20,0 to +250,0°C
- On-Off switch is light up green
- Setting of heat power from 1 to 99%
- Irreversible rapid shut-down when reaching the security temperature

Heating bath HB 1500-S

Equipment such as HB 1500, also with safety on separate microprocessor controller and an internal temperature sensor. Exceeding the safety temperature, the heating is switched off completely. The over-temperature shutdown >10K + reference deleted.



Code	Description	Packaging
LLG09234512	Heating bath HB 1500 / HB 1500-S, Type HB 1500, 2 litres, Min. operating temperature 20 RT + °C, Max. operating temperature 250 °C	1 pz
LLG09234510	Heating bath HB 1500 / HB 1500-S, Type HB 1500, 8 litres, Min. operating temperature 20 RT + °C, Max. operating temperature 250 °C	1 pz
LLG09234513	Heating bath HB 1500 / HB 1500-S, Type HB 1500-S, 8 litres, Min. operating temperature 20 RT + °C, Max. operating temperature 250 °C	1 pz
LLG09234511	Heating bath HB 1500 / HB 1500-S, Type HB 1500, 4 litres, Min. operating temperature 20 RT + °C, Max. operating temperature 250 °C	1 pz
LLG09234514	Heating bath HB 1500 / HB 1500-S, Type HB 1500-S, 4 litres, Min. operating temperature 20 RT + °C, Max. operating temperature 250 °C	1 pz
LLG09234515	Heating bath HB 1500 / HB 1500-S, Type HB 1500-S, 2 litres, Min. operating temperature 20 RT + °C, Max. operating temperature 250 °C	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

