

**Gas washing columns, FEP/PTFE****BOLA****Gas washing columns, FEP/PTFE**

PTFE top with two threaded connections. The standard PTFE frit has a pore size of approx. 3 μm and is screwed on the riser tube with an M8x1 thread. It can be exchanged with the PTFE gas distributor with fine bores which needs a lower primary pressure.

Temperature resistance: from -200°C to +250°C

-transparent

-unbreakable

-frit easily exchangeable



Code	Description	Packaging
LLG09110303	Gas washing columns, FEP/PTFE, Type Bottle, For external tube/pipe 8 mm, 1000 ml, Height 210 mm	
LLG09110308	Gas washing columns, FEP/PTFE, Type Column, For external tube/pipe 8 mm, 1000 ml, Height 700 mm	1 pz.
LLG09110301	Gas washing columns, FEP/PTFE, Type Bottle, For external tube/pipe 6 mm, 250 ml, Height 135 mm	
LLG09110302	Gas washing columns, FEP/PTFE, Type Bottle, For external tube/pipe 6 mm, 500 ml, Height 165 mm	
LLG09110306	Gas washing columns, FEP/PTFE, Type Column, For external tube/pipe 8 mm, 500 ml, Height 400 mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

