

**Magnetic stirrer MIX 1 XL****43953****Magnetic stirrer MIX 1 XL**

Wear-free 2mag motor drive, powerful, up to 40 liters, 70-2000 rpm, internal controller with 20 watts, incl. mixWATCH stirring bar monitoring with fully automatic trim of the max. Reachable speed for reliable and unattended operation, digital display, 4-step power setting for high power with regard to large and viscous media and reduced power for warming-free continuous operation, adjustable SoftStart for different volumes/viscosity, QuickStop against after-run and breakage, reliable sealed stainless-steel housing, 3 years warranty. Made in Germany.

Specifications

Operating conditions: -10°C to +50 °C at 95 % R.H. max.

Housing: Stainless steel

IP code: IP 64



Code	Description	Packaging
LLG06236109	Magnetic stirrer MIX 1 XL 100-2000 rpm, IP 64, SoftStart-function	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

