



Accessories for ThermoMixers and BlockThermostats

Hettich Benelux

Accessories for ThermoMixers and BlockThermostats



Code	Description	Packaging
LLG09777144	Accessories for ThermoMixers and BlockThermostats, Type BA 24, Anti-condensate plate for 1 block with micro tubes*	1 pz.
LLG09777145	Accessories for ThermoMixers and BlockThermostats, Type BA 96, Anti-condensate plate for 1 block with PCR and test plates**	1 pz.
LLG09777147	Accessories for ThermoMixers and BlockThermostats, Type BR 05, Assembly rack with stand for 0.5ml tubes	1 pz.
LLG09777146	Accessories for ThermoMixers and BlockThermostats, Type BH 10, Heating lid (200W), connects with cooling devices, for blocks for micro tubes and PCR plates	1 pz.
LLG09777142	Accessories for ThermoMixers and BlockThermostats, Type BT 01, Tempering cover for devices with 1 block, transparent, for high tubes	1 pz.
LLG09777143	Accessories for ThermoMixers and BlockThermostats, Type BT 02, Tempering cover for devices with 2 blocks, transparent, for high tubes	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

