



Safety suction systems AA, AC and AZ series

Safety suction systems AA, AC and AZ series

The DITABIS Aspiration Systems are made for the controlled aspiration of liquids from typical lab containers like cuvettes, tubes, microplates, petri dishes or roller bottles. They are also suited for the disposal of liquid biological or chemical waste and filtration of liquids by means of filter systems.

- Membrane vacuum pump
- 2L or 4L bottle incl. fast lock clutches, filter, ventilation insert and 3M silicone hose

Features AA series:

- Membrane vacuum pump turns on and off automatically
- AA 02/04 with 1 handle
- AA 24 with 4L bottle, 2 handles and Y connector
- Enlarge the automatic suction system by adding up to 4 working places and reduce your expenses for each user.

Features AC series:

- Membrane pump with manual switch-off
- 2L or 4L bottle

Features AZ series:

- Small membrane vacuum pump
- 2L or 4L bottle
- Bottle holder optional



Code	Description	Packaging
LLG09777220	Safety suction systems AA 02, automatic switch-on / switch-off, Bottle size 2000 ml, dim. (W x D x H) 160 x 230 x 470 mm, UV filters 17 / 100 nm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

