



LLG-Snap ring vials ND11, micro-vials, PP

LLG Labware

LLG-Snap ring vials ND11, micro-vials, PP

The LLG-Snap ring micro-vials are suitable for PFAS analyses.

- Made of PP
- Free from polyhaloolefins



| Code | Description | Packaging |
|-------------|---|-----------|
| LLG06267115 | LLG-Snap ring vials ND11, micro-vials, PP, Nominal capacity 0.3 ml, Dimensions (ØxH) 11.6 x 32 mm, Colour Transparent, Bottom shape Flat, Version Insert, 0.3 ml, Without labelling field | 100 pz |



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

