



Accessories for Thermo Scientific MaxQ HP Orbital Shakers

Accessories for Thermo Scientific MaxQ HP Orbital Shakers



Code	Description	Packaging
LLG04009429	Accessories for Thermo Scientific MaxQ HP Orbital Shakers, 125ml Flask Clip	1 pz.
LLG04009432	Accessories for Thermo Scientific MaxQ HP Orbital Shakers, 1l Flask Clip	1 pz.
LLG04009430	Accessories for Thermo Scientific MaxQ HP Orbital Shakers, 250/300ml Flask Clip	1 pz.
LLG04009427	Accessories for Thermo Scientific MaxQ HP Orbital Shakers, 25ml Flask Clip	1 pz.
LLG04009436	Accessories for Thermo Scientific MaxQ HP Orbital Shakers, 2.8l Flask Clip	1 pz.
LLG04009433	Accessories for Thermo Scientific MaxQ HP Orbital Shakers, 2l Flask Clip	1 pz.
LLG04009434	Accessories for Thermo Scientific MaxQ HP Orbital Shakers, 4l Flask Clip	1 pz.
LLG04009431	Accessories for Thermo Scientific MaxQ HP Orbital Shakers, 500ml Flask Clip	1 pz.
LLG04009428	Accessories for Thermo Scientific MaxQ HP Orbital Shakers, 50ml Flask Clip	1 pz.
LLG04009435	Accessories for Thermo Scientific MaxQ HP Orbital Shakers, 6l Flask Clip	1 pz.
LLG04009516	Accessories for Thermo Scientific MaxQ HP Orbital Shakers, Adjustable Angle Test Tube Holder with Rack for 10-13mm (3-7ml) Tubes	1 pz.
LLG04009517	Accessories for Thermo Scientific MaxQ HP Orbital Shakers, Adjustable Angle Test Tube Holder with Rack for 16-20mm (10-20ml) Tubes	1 pz.
LLG04009518	Accessories for Thermo Scientific MaxQ HP Orbital Shakers, Adjustable Angle Test Tube Holder with Rack for 21-25mm (30ml) Tubes	1 pz.
LLG04009519	Accessories for Thermo Scientific MaxQ HP Orbital Shakers, Adjustable Angle Test Tube Holder with Rack for 26-30mm (50ml) Tubes	1 pz.
LLG04009441	Accessories for Thermo Scientific MaxQ HP Orbital Shakers, Flask Clip Starter Kit with 25ml, 50ml, 250/300ml, 500ml, 1l and 2l Flask Clips (2 each)	1 pz.
LLG04009437	Accessories for Thermo Scientific MaxQ HP Orbital Shakers, Microplate Flask Clip	1 pz.
LLG04009521	Accessories for Thermo Scientific MaxQ HP Orbital Shakers, Microplate Rack for Deep Well Plates, 2 compartments	1 pz.



Code	Description	Packaging
LLG04009522	Accessories for Thermo Scientific MaxQ HP Orbital Shakers, Microplate Rack for Standard Plates, 3 compartments	1 pz.
LLG04009532	Accessories for Thermo Scientific MaxQ HP Orbital Shakers, Test Tube Rack for 10-13mm (3-7ml) Tubes	1 pz.
LLG04009533	Accessories for Thermo Scientific MaxQ HP Orbital Shakers, Test Tube Rack for 16-20mm (10-20ml) Tubes	1 pz.
LLG04009514	Accessories for Thermo Scientific MaxQ HP Orbital Shakers, Test Tube Rack for 21-25mm (30ml) Tubes	1 pz.
LLG04009515	Accessories for Thermo Scientific MaxQ HP Orbital Shakers, Test Tube Rack for 26-30mm (50ml) Tubes	1 pz.
LLG04009353	Accessories for Thermo Scientific MaxQ HP Orbital Shakers, Tube Inserts, perforated to hold 3x50ml, 3x20ml and 8x 1.5/2.0ml Tubes	5 pz.
LLG04009438	Accessories for Thermo Scientific MaxQ HP Orbital Shakers, Utility Tray with Rubber Pad, 229x457mm	1 pz.
LLG04009439	Accessories for Thermo Scientific MaxQ HP Orbital Shakers, Utility Tray with Rubber Pad, 457x457mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

