

**LLG pH meter 7****LLG Labware****LLG pH meter 7**

The LLG pH Meter 7 is an easy-to-use, robust table-top instrument for measuring pH and redox potential. The large color-touch-display and the intuitive, language-independent user interface ensure fast and easy use. Particularly suitable for routine applications in laboratories, production plants and schools.

-Available as set with gel combination electrode ("Made in Germany") with temperature sensor

- Large color-touch-display with language-independent graphical user interface
- Automatic temperature compensation
- Memory for 500 values
- 3-point calibration
- Reminder function for upcoming calibrations
- BNC connection for the electrode
- RS232 interface for data transfer
- Electrode holder can be mounted left and right

Scope of supply:

LLG pH meter 7: pH meter, electrode holder, protective cover, power supply unit

LLG-pH Meter 7 set: pH meter, combination electrode, electrode holder, protective cover, power supply unit

Specifications

pH range: 0 ... 14

pH resolution: 0.01

pH accuracy: $\pm 0,01$

Temperature range electrode: 0 ... 80 °C

Dimensions (W x D x H): 180 x 240 x 80 mm

Weight: 0.7 kg

Power supply: 9 V DC/1A

IP code: IP 54

Warranty: 3 years



Code	Description	Packaging
LLG06263600	LLG-pH Meter 7	1 pz.
LLG06263601	set LLG-pH Meter 7	1 pz.
LLG06263607	set LLG-pH Meter 7 (UK plug)	1 pc



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

