



LLG-Storage Boxes, PP

LLG Labware

LLG-Storage Boxes, PP



Code	Description	Packaging
LLG04001528	LLG-Storage Boxes, PP, For 5/10/20 ml, Max. no. of sample positions 25, Dimensions (WxDxH) 130 x 130 x 102 mm, Colour Blue, Material PP, Fits for ND20	1 pz
LLG09405760	LLG-Storage Boxes, PP, For 1.5/2.0 ml, Max. no. of sample positions 16, Dimensions (WxDxH) 67 x 67 x 45 mm, Colour Green, Fits for ND8, 9, 10, 11	1 pz
LLG09405761	LLG-Storage Boxes, PP, For 1.5/2.0 ml, Max. no. of sample positions 16, Dimensions (WxDxH) 67 x 67 x 45 mm, Colour Pink, Fits for ND8, 9, 10, 11	1 pz
LLG09405765	LLG-Storage Boxes, PP, For 1.5/2.0 ml, Max. no. of sample positions 16, Dimensions (WxDxH) 67 x 67 x 45 mm, Colour Transparent, Fits for ND8, 9, 10, 11	1 pz



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

