

**Floating racks PrepSafe for microcentrifuge tubes, PP****Bel-Art Products****Floating racks PrepSafe for microcentrifuge tubes, PP**

Unique tabs protect your samples from loss and contamination. Place tubes in slots and slide under tabs to keep tubes in place and caps closed when under pressure, floating or submerged. Use in water baths, ice baths or on the benchtop period.

- Holds up to twelve 1.5ml or 2.0ml microcentrifuge tubes in numbered locations
- Raised handle for easy lifting without agitating samples; stable feet keep racks upright on the bench
- 95mm outer diameter fits most 1000ml beakers
- With removable vortexing attachment, this allows samples to be vortexed simultaneously
- Steam autoclavable at 121°C



Code	Description	Packaging
LLG09400995	Floating racks PrepSafe for microcentrifuge tubes, PP, Colour Natural	1 pz.
LLG09400998	Floating racks PrepSafe for microcentrifuge tubes, PP, Colour Neon Yellow	1 pz.
LLG09400997	Floating racks PrepSafe for microcentrifuge tubes, PP, Colour Pink	1 pz.
LLG09400996	Floating racks PrepSafe for microcentrifuge tubes, PP, Colour Turquoise	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

