

**LLG pH electrode****LLG Labware****LLG pH electrode**

Low maintenance combination electrode with integrated temperature sensor.
Robust design with plastic shaft, protective basket and fixed cable connection.
Made in Germany.

Specifications**Temperature range:****Temperature sensor: 0 ... 80 °C****NTC 30 kW****pH range: 0 ... 14****Diaphragm: 2 x ceramic****Membrane shape: Cylindrical****Electrolyte: Gel, 3M****Min. immersion depth: 15 mm****Cable length: 1 m****Length/Ø: 120/12 mm****Electrode connection:****Temperature sensor connection: BNC****Cinch 3 mm**

Code	Description	Packaging
LLG06263602	pH-electrode for LLG-pH meter 7	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

