

**Cryogenic boxes, slip lid with adaptable height, PP, 133 x 133****Ratiolab****Cryogenic boxes, slip lid with adaptable height, PP, 133 x 133**

- Height adaptable from 50 mm to 75 mm, by turning the lid 90°
- 9 x 9 array, 10 x 10 array
- Numerical coded places (except black box)
- Drain bores at the bottom for dew liquids
- Autoclavable, temperature resistant from -90 °C to 121 °C



Code	Description	Packaging
LLG09405847	Cryogenic boxes, slip lid with adaptable height, PP, 133 x 133, For tubes diam. 11.5 mm, Colour Natural, Array 10 x 10	5 pz.
LLG09405848	Cryogenic boxes, slip lid with adaptable height, PP, 133 x 133, For tubes diam. 11.5 mm, Colour Yellow, Array 10 x 10	5 pz.
LLG09405849	Cryogenic boxes, slip lid with adaptable height, PP, 133 x 133, For tubes diam. 11.5 mm, Colour Red, Array 10 x 10	5 pz.
LLG09405850	Cryogenic boxes, slip lid with adaptable height, PP, 133 x 133, For tubes diam. 11.5 mm, Colour Green, Array 10 x 10	5 pz.
LLG09405851	Cryogenic boxes, slip lid with adaptable height, PP, 133 x 133, For tubes diam. 11.5 mm, Colour Blue, Array 10 x 10	5 pz.
LLG09405852	Cryogenic boxes, slip lid with adaptable height, PP, 133 x 133, For tubes diam. 11.5 mm, Colour Black, Array 10 x 10	5 pz.
LLG09405855	Cryogenic boxes, slip lid with adaptable height, PP, 133 x 133, For tubes diam. 11.5 mm, Colour Natural, yellow, red, green, blue, Array 10 x 10	5 pz.
LLG09405840	Cryogenic boxes, slip lid with adaptable height, PP, 133 x 133, For tubes diam. 13 mm, Colour Natural, Array 9 x 9	5 pz.
LLG09405841	Cryogenic boxes, slip lid with adaptable height, PP, 133 x 133, For tubes diam. 13 mm, Colour Natural, yellow, red, green, blue, Array 9 x 9	5 pz.
LLG09405842	Cryogenic boxes, slip lid with adaptable height, PP, 133 x 133, For tubes diam. 13 mm, Colour Blue, Array 9 x 9	5 pz.
LLG09405843	Cryogenic boxes, slip lid with adaptable height, PP, 133 x 133, For tubes diam. 13 mm, Colour Yellow, Array 9 x 9	5 pz.
LLG09405844	Cryogenic boxes, slip lid with adaptable height, PP, 133 x 133, For tubes diam. 13 mm, Colour Green, Array 9 x 9	5 pz.
LLG09405845	Cryogenic boxes, slip lid with adaptable height, PP, 133 x 133, For tubes diam. 13 mm, Colour Red, Array 9 x 9	5 pz.
LLG09405846	Cryogenic boxes, slip lid with adaptable height, PP, 133 x 133, For tubes diam. 13 mm, Colour Black, Array 9 x 9	5 pz.





CARLO ERBA Reagents S.A.S.
Chaussée du Vexin – 27106 Val
de Reuil cedex
N° TVA: FR 63391048824
Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

