



Vortexer LLG-uniTEXER with universal attachment

Vortexer LLG-uniTEXER with universal attachment One head - multiple uses

Vortexer with single carrier head, accommodating the most popular consumables and tubes:
for 1 x assay plate, 2 x 50 ml conical tubes, 2 x 15 ml conical tubes, 2 x 5 ml conical tubes, 4 x 1.5/2 ml microcentrifuge tubes, 6 x 0.5 ml and 24 x 0.2 ml microcentrifuge tubes.

- Touch function for short-term operation
- Very easy to clean
- Elastomeric feet ensure excellent stability on the bench
- Modern design
- Robust housing made of ABS

Scope of supply: vortexer, 1 x universal attachment

Specifications

Type of movement: orbital/vortex

Orbital diameter: 3.7 mm

Speed: 1000, 2000 und 3000 min⁻¹

Dimensions (W x D x H): 173 x 198 x 198 mm

Weight: 3.8 kg

Power supply: 220 V, 50/60 Hz

IP code: IP 20

Warranty: 3 years



Code	Description	Packaging
LLG06261826	uniTEXER, vortexer with universal attachment, with EU plug	1 pezzo



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

