



Hand held thermometer Type 13100

Hand held thermometer Type 13100

Digital temperature measuring device with 2 input channels, incl. 1 sensor with flat plug (measuring range up to approx. 150 °C, cable length: 3 m). For measuring sensors type K (NiCr-Ni). Measuring range: -200 ...1370 °C / -328 ... 2498 °F. Resolution: 0.1 °C/F. Accuracy: ± 0.1 %. With hold function, with max/min memory and automatic switch-off after 30 min. With PC interface and bag. Battery: 1 x 9 V block. Dimensions: 184 x 64 x 30 mm



Code	Description	Packaging
LLG09236719	Hand held thermometer Type 13100, 90° Surface sensor with handle type K	1 pz.
LLG09236715	Hand held thermometer Type 13100, Hand held thermometer Type 13100	1 pezzo
LLG09236717	Hand held thermometer Type 13100, Insertion sensor with handle type K	1 pz.
LLG09236716	Hand held thermometer Type 13100, Power adapter (100 ... 240 V, 50 ... 60 Hz)	1 pz.
LLG09236718	Hand held thermometer Type 13100, Surface sensor with handle type K	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

