

**LLG-Electrical pipette Plus****LLG Labware****LLG-Electrical pipette *Plus***

Ergonomic and precise - the cordless and lightweight LLG-Electrical pipette *Plus* enables fast and exact pipetting with glass and plastic pipettes.

- UV resistant housing
- Fatigue-proof pipetting due to ergonomically shaped handle
- Suitable for left and right handers
- Exchangeable filters protect the device mechanism from contamination
- LED display with speed display: High, Low or Gravity mode (for free liquid flow without drive operation)
- 5 different coloured autoclavable nose pieces and labelling area enables customisation of the pipette
- Extremely quiet and powerful drive
- Low-Battery LED indicator

Scope of supply: Charger (5 V micro USB), Li-Ion battery, 5 coloured nose pieces, 3 membrane filters, adapter for 1 ml pipettes, benchtop rack, magnetic wall mounting

Specifications

Membrane filters: PTFE/Ø 25 mm /0.45 µm

Capacity: 1 to 100 ml, glass and plastic pipettes

Pipette attachment: Polypropylene, autoclavable

Dimensions: 15 x 152 x 41 mm

Weight: approx. 200 g

Power supply: rechargeable Li-Ion battery, 8 h continuous operation

Warranty: 1 year



Code	Description	Packaging
LLG09281041	Electrical pipettePlus, Electrical pipette Plus	1 pz
LLG09281005	Electrical pipettePlus, Membrane filter, 0.45µm, PTFE, non-sterile	5 pz
LLG06267647	Electrical pipettePlus, Replacement accu	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

