

**Halogen Moisture Analyser DBS 60-3****Kern & Sohn****Halogen Moisture Analyser DBS 60-3**

Modern moisture analyser with 10 memories and graphic display

- Halogen quartz glass heater 400 W
- Backlit LCD display, digit height 15 mm: Previous drying time, current temperature, unit for displaying the results, current moisture content in %, drying mode/Status display drying
- Observation window over the sample, useful during initial adjustment
- Internal memory for automatic sequence of 10 drying programs and 100 drying processes carried out
- The last value measured remains on the display until it is replaced by a new measurement
- Sample description for up to 99 samples, 2 digits, freely programmable, and is printed in the measuring protocol
- Date and time display as standard
- 10 sample plates included
- Table of applications: There are many practical examples in the operating instructions

**Specifications**

**Readability: 0.01 % or 0.001 g**

**Weighing range, max: 60 g**

**Reproducibility, weight of sample 2 g: 0.15 %**

**Reproducibility, weight of sample 10 g: 0.02 %**

**Temperature range: +50 to +200 °C**

**Dimensions (W x D x H): 204 x 336 x 167 mm**



| Code        | Description   | Packaging |
|-------------|---|-----------|
| LLG09901441 | Halogen Moisture Analyser DBS 60-3, Type DBS 60-3, max. load 60 g | 1 pz      |



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

