

**Halogen Moisture Analyser DBS 60-3****Kern & Sohn****Halogen Moisture Analyser DBS 60-3**

Modern moisture analyser with 10 memories and graphic display

- Halogen quartz glass heater 400 W
- Backlit LCD display, digit height 15 mm: Previous drying time, current temperature, unit for displaying the results, current moisture content in %, drying mode/Status display drying
- Observation window over the sample, useful during initial adjustment
- Internal memory for automatic sequence of 10 drying programs and 100 drying processes carried out
- The last value measured remains on the display until it is replaced by a new measurement
- Sample description for up to 99 samples, 2 digits, freely programmable, and is printed in the measuring protocol
- Date and time display as standard
- 10 sample plates included
- Table of applications: There are many practical examples in the operating instructions

Specifications

Readability: 0.01 % or 0.001 g

Weighing range, max: 60 g

Reproducibility, weight of sample 2 g: 0.15 %

Reproducibility, weight of sample 10 g: 0.02 %

Temperature range: +50 to +200 °C

Dimensions (W x D x H): 204 x 336 x 167 mm



Code	Description	Packaging
LLG09901441	Halogen Moisture Analyser DBS 60-3, Type DBS 60-3, max. load 60 g	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

